

# **MEDICAL LCD MONITOR**

## **USER'S GUIDE**

---

Before connecting, operating or adjusting this product, please read this instruction booklet carefully and completely.

# MEDICAL LCD MONITOR

FS-L190\*D(19")

FS-L190\*DT(19")

FS-L240\*D(24")

FS-L240\*DT(24")

FS-L260\*D(26")

FS-L320\*D(32")

FS-L420\*D(42")

## Model definition

F S- L 190 \* D T

1. FS : Monitor manufacturer.
2. L : Panel manufacturer.
3. 190: Display size.
4. \* : Signal input option.
5. D : Medical grade.
6. T : Touch screen installed.

# Table of Contents

Symbol Definitions .....	2
Safety Instructions .....	4
Caution .....	7
FCC Information .....	10
Parts .....	15
● FS-L190*D / FS-L190*DT .....	15
● FS-L240*D / FS-L240*DT .....	16
● FS-L260*D .....	17
● FS-L320*D .....	18
● FS-L420*D .....	19
Connector .....	20
● FS-L190*D / FS-L240*D / FS-L260*D / FS-L320*D .....	20
● FS-L190*DT / FS-L240*DT .....	21
● FS-L420*D .....	22
Mechanical Product Drawing .....	23
● FS-L190*D / FS-L190*DT .....	23
● FS-L240*D / FS-L240*DT .....	24
● FS-L260*D .....	25
● FS-L320*D .....	26
● FS-L420*D .....	27
Control .....	28
Power management .....	32
OSD .....	33
Standard Signal table .....	48
Signal connector Pin Assignments .....	49
Specification .....	52
● FS-L190*D / FS-L190*DT .....	52
● FS-L240*D / FS-L240*DT .....	53
● FS-L260*D .....	54
● FS-L320*D .....	55
● FS-L420*D .....	56

## Symbol Definitions

The following symbols appear on the product, its labels, or the product package. Each symbol carries a special definition, as defined below.



Dangerous : High Voltage.



Consult accompanying documents.



Direct Current.



Indicates protective earth ground.



DC Power control switch.

**SN**

Serial Number.



Top-Bottom.



Fragile.



Do not get wet.

**3**

Maximum Stacking.(19"/24"/26"/32")

**2**

Maximum Stacking.(42")

**CE**

Indicates proof of conformity to applicable European Economic Community Council directives and to harmonized standards published in the official journal of the European Communities.



Medical LCD monitor is in accordance with UL 60601-1 and CAN/CSA C22.2 No.601.1 in regards to electric shock, fire hazards, and mechanical hazard.

**FCC**

Tested to comply with FCC Class B standard.



This symbol indicates that the waste of medical LCD monitor must not be disposed as unsorted municipal waste and must be collected separately. Please contact the manufacturer or other authorized disposal company to decommission your medical LCD monitor.

## Safety Instructions

### On Safety

1. Before connecting the AC power cord to the DC adapter outlet make sure the voltage designation of the DC adapter corresponds to the local electrical supply.
2. Never insert anything metallic into the cabinet openings of the medical LCD monitor. Doing so may create the danger of electric shock.
3. To reduce the risk of electric shock, do not remove cover.  
No user-serviceable parts inside. Only a qualified technician should open the case of the medical LCD monitor.
4. Never use your medical LCD monitor if the power cord has been damaged.  
Do not allow anything to rest on the power cord, and keep the cord away from areas where people can trip over it.
5. Be sure to hold the plug, not the cord, when disconnecting the medical LCD monitor power cord from an electric socket.
6. Unplug your medical LCD monitor power cord when it is going to be left unused for an extended period of time.
7. Unplug your medical LCD monitor power cord from the AC outlet before any service.
8. If your medical LCD monitor does not operate normally, in particular, if there are any unusual sounds or smells coming from it, unplug it immediately and contact an authorized dealer or service center.
9. Please contact the manufacturer if the set should be installed in an inaccessible area.

### Warning

Do not touch input or output connectors and the patient simultaneously.

## **Warning**

This medical LCD monitor is intended for connection to input/output signals and other connectors that comply with relevant IEC standard (e.g., IEC60950 for IT equipment and IEC60601 series for medical electrical equipment).

In addition, all such combination-system shall comply with the standard IEC 60601-1-1, safety requirements for medical electrical systems. Any person who has formed a combination-system is responsible for the system to comply with the requirements of IEC 60601-1-1.

If in doubt, contact a qualified technician or your local representative.

## **On installation**

1. Openings in the medical LCD monitor cabinet are provided for ventilation.  
To prevent overheating, these openings should not be blocked or covered. If you put the medical LCD monitor in a bookcase or some other enclosed space, be sure to provide adequate ventilation.
2. Put your medical LCD monitor in a location with low humidity and a minimum of dust.
3. Do not expose the medical LCD monitor to rain or use it near water (in kitchens, near swimming pools, etc.). If the medical LCD monitor accidentally gets wet, unplug it and contact an authorized dealer immediately. You can clean the medical LCD monitor with a damp cloth if necessary, but be sure to unplug the medical LCD monitor first.
4. Place your medical LCD monitor near an easily accessible AC outlet.
5. High temperature can cause problems. Don't use your medical LCD monitor in direct sunlight and keep it away from heaters, stoves, fireplaces, and sources of heat.
6. Don't place your medical LCD Monitor on an unstable stand, Medical LCD monitor may malfunction or fall.
7. This medical LCD monitor should not topple over when tilted at a 5 degree angle, in any position, during NORMAL USE, excluding transport.
8. In the position specified for transport, Medical LCD monitor shall not overbalance when tilted at a 10 degree angle.

#### **Environmental Conditions for operation and Storage**

- Temperature range within 0°C to 40°C(operation), -20°C to 60°C(storage)
  - Relative humidity range 10% to 85%
- Atmospheric pressure range within 500 to 1060hPa.

#### **Intended Use**

- This Medical LCD Monitor is an accessory intended for use with Medical Equipment to display alphabetical, numerical and graphical data.

## CAUTION



This symbol alerts the user that important literature concerning the operation of this unit has been included. Therefore, it should be read carefully in order to avoid potential problems.



This symbol warns user that un-insulated voltage within the unit may have sufficient magnitude to cause electrical shock. Therefore, it is dangerous to make contact with any part inside the unit. To reduce the risk of electric shock, DO NOT remove cover (or back).

There are no user serviceable parts inside. Refer servicing to qualified service personnel.

To prevent fire or shock hazards, do not expose this unit to rain or moisture. Also, do not use this unit's polarized plug with an extension cord receptacle or other outlets unless the prongs can be fully inserted. The display is designed to meet the medical safety requirements for a patient vicinity device.

This device may not be used in connection with life support equipment.



Underwriters Laboratories (UL) Classification:

UL safety Compliance:

This medical LCD monitor is U.L. Classified WITH RESPECT TO ELECTRIC SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH UL 60601-1/CAN/CSA C22.2 NO. 601.1



### **EEC Safety Compliance**

This medical LCD monitor unit meets the requirements of EN-60601-1 so as to conform to the Medical Device Directive 93/42/EEC (general safety information).

Use 120V rating 5-15P type plug only in the U.S

This medical LCD monitor complies to the above standards only when used with the supplied medical grade power supply.

19"(FS-L190\*D / FS-L190\*DT)

JMW190KA1200F02(BRIDGE POWER CORP.)

24"/26"(FS-L240\*D / FS-L240\*DT / FS-L260\*D)

JMW1150KA2400F04(BRIDGE POWER CORP.)

32"(FS-L320\*D) JMW1180KA2400F01(BRIDGE POWER CORP.)

Caution: Make sure the power cord is the correct type that is required in your area.

This medical LCD monitor has a universal power supply that allows operation in either 100-120V AC or 200-240V AC voltage areas (no user adjustment is required).

Use the proper power cord with correct attachment plug type. If the power source is 120 V AC, use a power cord which is a Hospital Grade Power Cord with NEMA 5-15 style plug, labeled for 125 volts AC with UL and C-UL approvals. If the power source is a 240 V AC supply, use the tandem (T blade) type attachment plug with ground conductor power cord that meets the respective European country's safety regulations.

The hospital-grade plug for medical products intended for use in Denmark has DEMKO approval and is rated 13 amps at 250Vac. Plug is recommended for use in medical applications and specifications are being added to the standard SB 107-2-D1. Plug mates with maker's Danish hospital-grade socket. Hospital sockets have slightly different shaped openings allowing only the hospital plug, not the standard Danish plug, to be inserted, to protect the ac circuit in specific medical settings.

### **Recycling**



Follow local governing ordinances and recycling plans regarding the recycling or disposal of this equipment.

## **Cleaning Instructions**



Follow your hospital protocol for the handling of blood and body fluids. Clean the display with a diluted mixture of mild detergent and water. Use a soft towel or swab.

Use of certain detergents may cause degradation to the labels and plastic components of the product.

Consult cleanser manufacturer to see if agent is compatible with it.

Do not allow liquid enter the display.

## **Servicing**

Do not attempt to service the medical LCD monitor yourself, as opening or removing covers may expose you to dangerous voltages or other hazards, and will void the warranty. Refer all servicing to qualified service personnel.

Unplug the medical LCD monitor from its power source and refer servicing to qualified personnel under the following conditions:

- If the power cord or plug is damaged or frayed.
- If liquid has been spilled into the medical LCD monitor.
- If objects have fallen into the medical LCD monitor.
- If the medical LCD monitor has been exposed to rain or moisture.
- If the medical LCD monitor has been subjected to excessive shock by being dropped.
- If the cabinet has been damaged.
- If the medical LCD monitor seems to be overheated.
- If the medical LCD monitor emits smoke or abnormal odor.
- If the medical LCD monitor fails to operate in accordance with the operating instructions.

## **Accessories**

Use only accessories specified by the manufacturer, or sold with the medical LCD monitor.

## **Classification**

- Protection against electric shock : Class I including AC/DC Adapter
- Applied Parts : No Applied Parts
- Degree of safety in the presence of flammable anesthetics mixture with air or with oxygen or with nitrous oxide.  
Not suitable for use in the presence of a flammable anesthetics mixture with oxygen or with nitrous oxide.
- Mode of operation : Continuous.

## FCC Information

This medical LCD monitor unit has been tested and found to comply with the limits of a Class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against interference. This monitor can radiate radio frequency energy and, if not installed and used in accordance with the instructions, it may interfere with other radio communications equipment. There is no guarantee that interference will not occur in a particular installation.

If this equipment is found to cause harmful interference to radio or television reception, the user is encouraged to try to correct the interference by carrying out one or more of the following measures:

1. Reorient or relocate the receiving antenna.
2. Increase the distance between the medical LCD monitor and the subject of interference.
3. Plug the monitor into an outlet on a different electrical circuit than that to which the subject of interference is connected.
4. Consult the dealer or an experienced radio/TV technician for help.

### NOTICES TO USER

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

### FCC WARNING

This medical LCD monitor generates or uses radio frequency energy. Changes or modifications to this medical LCD monitor may cause harmful interference unless the modifications are expressly approved in the instruction manual. The user could lose authority to operate this equipment if an unauthorized change or modification is made.

## 1. Guidance and manufacturer's declaration - electromagnetic emissions

The medical LCD monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the medical LCD monitor should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment -guidance
RF Emissions CISPR 11	Group 1	The medical LCD monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF Emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	D	The medical LCD monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Voltage fluctuations IEC 61000-3-3	Complies	

## 2. Guidance and manufacturer's declaration - electromagnetic immunity

This medical LCD monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the medical LCD monitor should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment -guidance
Electrostatic discharge(ESD) IEC61000-4-2	6kV Contact 8kV air	6kV Contact 8kV air	Floors should be wood,concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	2kV for power supply lines 1kV for input/output lines	2kV for power supply lines 1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.

Surge IEC 61000-4-5	1kV differential mode 2kV common mode	1kV differential mode 2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz	3 Vrms 150 kHz to 80MHz	<p>Portable and mobile RF communications equipment should be used no closer to any part of the medical LCD monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance : <math>d</math></b></p> $d = \left[ \frac{3.5}{V_f} \right] \sqrt{P}$ <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W)</p>

### 3. Guidance and manufacturer's declaration - electromagnetic immunity

<p>This medical LCD monitor is intended for use in the electromagnetic environment specified below. The customer or the user of monitor should assure that it is used in such an environment.</p>			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment -guidance
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3.0A/m	3.0A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<p>&lt;5 % <math>UT</math> (&gt;95 % dip in <math>UT</math>) for 0.5 cycle</p> <p>40 % <math>UT</math> (60 % dip in <math>UT</math>) for 5 cycle</p> <p>70 % <math>UT</math> (30 % dip in <math>UT</math>) for 25 cycle</p> <p>&lt;5 % <math>UT</math> (&lt;95 % dip in <math>UT</math>) for 5 sec.</p>	<p>&lt;5 % <math>UT</math> (&gt;95 % dip in <math>UT</math>) for 0.5 cycle</p> <p>40 % <math>UT</math> (60 % dip in <math>UT</math>) for 5 cycle</p> <p>70 % <math>UT</math> (30 % dip in <math>UT</math>) for 25 cycle</p> <p>&lt;5 % <math>UT</math> (&lt;95 % dip in <math>UT</math>) for 5 sec.</p>	<p>Main power quality should be that of a typical commercial or hospital environment. If the user of monitor requires continued operation during power mains interruptions, it is recommended that monitor be powered from an uninterruptible power supply or a battery.</p> <p>NOTE : <math>UT</math> the A.C. mains voltage prior to application of the test level.</p>

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment -guidance
Radiated RF IEC 61000-4-3	3 V/m 80.0 MHz to 2.5 GHz	3 V/m 80.0 MHz to 2.5 GHz	<p><b>Recommended separation distance</b></p> <p>80MHz to 800MHz</p> $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ <p>80MHz to 2.5GHz</p> $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, It should be less than the compliance level in each frequency range.</p>

#### 4. Recommended separation distances between portable and mobile RF communications equipment and this medical LCD monitor.

- The medical LCD monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled.
- The customer or the user of the monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the medical LCD monitor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance according to frequency of transmitter[m]		
	150kHz to 80MHz	80MHz to 800MHz	800MHz to 2.5GHz
	$d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$	$d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	$d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
	V1=3Vrms	E1=3V/m	E1=3V/m

0.01	0.116	0.116	0.2333
0.1	0.368	0.3687	0.7378
1	1.166	1.1660	0.2333
10	3.687	3.6872	0.7375
100	11.660	11.6600	23.333


For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $p$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1) At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Parts


### 19" FS-L190\*D / FS-L190\*DT Monitor




### Accessories




User Manual




D-SUB Cable




AC-DC Adaptor  
(JMW190KA1200F02)



S-Video (Y/C) Cable  
(Option)




DVI Cable




BNC Cable



AC Power cord  
(Hospital Grade)



SCREW FH M3X6



SCREW BH M4X10


▶ 24" FS-L240\*D / FS-L240\*DT Monitor




▶ Accessories




User Manual




D-SUB Cable



AC-DC Adaptor  
(JMW1150KA2400F04)



S-Video (Y/C) Cable  
(Option)




DVI Cable



BNC Cable




AC Power cord  
(Hospital Grade)



DC Cable Terminal  
Male / female (Option)




SCREW FH M3X6



SCREW BH M4X10


● 26" FS-L260\*D Monitor




● Accessories




User Manual




D-SUB Cable




AC-DC Adaptor  
(JMW1150KA2400F04)



S-Video (Y/C) Cable  
(Option)




DVI Cable




BNC Cable




AC Power cord  
(Hospital Grade)



DC Cable Terminal  
Male / female (Option)



SCREW FH M3X6



SCREW BH M4X10


● 32" FS-L320\*D Monitor




● Accessories




User Manual




D-SUB Cable




AC-DC Adaptor  
(JMW1180KA2400F01)



S-Video (Y/C) Cable  
(Option)




DVI Cable




BNC Cable




AC Power cord  
(Hospital Grade)



DC Cable Terminal  
Male / female (Option)




SCREW FH M3X6



SCREW BH M4X10



Remote controller



Battery(AAA)


▶ 42" FS-L420\*D Monitor



▶ Accessories




User Manual




D-SUB Cable




S-Video (Y/C) Cable  
(Option)




DVI Cable



BNC Cable



AC Power cord  
(Hospital Grade)




Remote controller




Battery(AAA)

## Connector


► 19" FS-L190\*D / FS-L240\*D / FS-L260\*D / FS-L320\*D Input connector




► Dual DVI input




▶ FS-L190\*DT / FS-L240\*DT input connector




▶ Dual DVI input



▶ FS-L420\*D Input connector




▶ Dual DVI input




## Mechanical Product Drawing

### FS-L190\*D / FS-L190\*DT Dimension




Front view



Rear view (Without I/O Cover)




Rear view (Installed I/O Cover)





Unit : mm

▶ FS-L240\*D / FS-L240\*DT Dimension



Front view



Rear view (Without I/O Cover)




Rear view (Installed I/O Cover)




Unit : mm

FS-L260\*D Dimension


### Front view



### Rear view (Without I/O Cover)




### Rear view (Installed I/O Cover)




Unit : mm

▶ FS-L320\*D Dimension


Front view



Rear view (Without I/O Cover)




Rear view (Installed I/O Cover)





Unit : mm

▶ FS-L420\*D Dimension

Front view




Rear view



Unit : mm

## Control

### OSD Button



An 8 button keypad, located in bottom right corner on the front of the display, allows the user to make adjustments to various display parameters using the On Screen Display (OSD) system.

### ● Power Indicating LED

Normal mode : Green

Standby mode : Amber

Monitor Off : Off

Note : The LED state of Normal Mode and Monitor Off may be reversed.

This can be change upon request.



### ● On-Screen Display (OSD) Function Button

1. POWER : Soft power turns the monitor ON or OFF.
2. MENU : Used to activate to the OSD menu and exit from main menu or sub menu.
3. PIP : Enables PIP (picture in Picture) function. Select PIP, PBP1 or PBP2.
4. UP ( $\blacktriangle$ ) : With the OSD deactivated, it is a Hot Key for increasing brightness.  
With the OSD activated, moves the cursor upward.
5. DOWN ( $\blacktriangledown$ ) : With the OSD deactivated, it is a Hot Key for decreasing brightness.  
With the OSD activated, moves the cursor downward.

- 6. PLUS (+) : With the OSD deactivated, it is a Hot Key for increasing contrast.  
With the OSD activated, enter sub menu and increases the adjustment of the selected function.
- 7. MINUS (-) : With the OSD deactivated, it is a Hot Key for decreasing contrast.  
With the OSD activated, it decreases the adjustment of the selected function.
- 8. INPUT : With the OSD deactivated, by pressing down for over 1 sec., it is a Hot Key for auto-adjustment control under DSUB ANALOG / RGBs signals.  
With the OSD activated, it changes the displayed signal source.

▶ Remote button function

Note: Remote control is available for the FS-L320\*D and FS-L420\*D models only.




1. SOURCE : Change the display signal source.
2. POWER : Soft power turns the monitor ON or OFF.
3. UP (▲) : With the OSD deactivated, it is a Hot key for increasing brightness.  
With the OSD activated, moves the cursor upward.
4. DOWN (▼) : With the OSD deactivated, it is a Hot Key for decreasing brightness.  
With the OSD activated, moves the cursor downward.
5. MINUS (-) : With OSD deactivated, it is a Hot Key for decreasing contrast.  
With the OSD activated, it decreases the adjustment of the selected function.
6. PLUS (+) : With the OSD deactivated, it is Hot Key for increasing contrast.  
With the OSD activated, it decreases the adjustment of the selected function.
7. MENU : Used to activate to OSD menu and exit from main menu or sub menu.
8. AUTO : Fit to the most appropriate screen on the D-SUB Analog signal.
9. PIP : Enables PIP(picture in Picture) function.  
Select PIP,PBP1,PBP2
10. MUTE : Sound muted.
11. SWAP : Swaps the position of the Primary and Secondary images.

## ▶ GPIO


### 1. Primary and Secondary Swap

Pressing the swap button swaps the primary and secondary image.




### 2. PIP,PBP1,PBP2 Select

The selection of secondary image position button pressed.




### 3. Record Indicator

The record indicator is displayed while a contact closure to the ground pin is present. The indicator is removed when the contacts are opened. The record indicator is displayed in the monitor's top left corner.



### 4. GND

Common ground



## Power management

This monitor does not adhere to the VESA DPMS standard when no signal is present on the video input.


Model	Status	LED sign	Power Consumption
FS-L190*D	Normal mode	Green on	<60W
	Standby mode	Green blinking	<20W
FS-L190*DT	Normal mode	Green on	<60W
	Standby mode	Green blinking	<20W
FS-L240*D	Normal mode	Green on	<100W
	Standby mode	Green blinking	<20W
FS-L240*DT	Normal mode	Green on	<100W
	Standby mode	Green blinking	<20W
FS-L260*D	Normal mode	Green on	<130W
	Standby mode	Green blinking	<20W
FS-L320*D	Normal mode	Green on	<150W
	Standby mode	Green blinking	<20W
FS-L420*D	Normal mode	Green on	<260W
	Standby mode	Green blinking	<20W

Note : The LED state of Normal Mode and Monitor Off may be reversed.  
This can be changed upon request.

# OSD




## ADJUST




1. BRIGHTNESS  
Increase or decrease the brightness. (Range : 0~100)
2. CONTRAST  
Increase or decrease the contrast. (Range : 0~100)
3. CLOCK  
Increase or decrease the sampling frequency. (Range : 0~100)
4. PHASE  
Increase or decrease the Phase level. (Range : 0~100)
5. BACKLIGHT  
Adjust backlight dimming level. (Range : 0~100)
6. AUTO ADJUST  
Fit to the most appropriate screen on the D-SUB Analog / RGBs signal.

## COLOR TEMP




1. MODE  
Change the color mode C1(Reddish,6500K), C2(Bluish,9300K), USER(7200K)
2. RED  
Red balance. (Only works with USER Mode) (Range : 0~100)
3. GREEN  
Green balance. (Only works with USER Mode) (Range : 0~100)
4. BLUE  
Blue balance. (Only works with USER Mode) (Range : 0~100)

## IMAGE




1. IMAGE SIZE  
Change the image size. (Full, Fill aspect, 1:1, Normal)
2. H POSITION  
Adjust the horizontal position of the displayed source image.(Range : 0~100)
3. V POSITION  
Adjust the vertical position of the displayed source image. (Range : 0~100)
4. GAMMA  
Adjust GAMMA value (VIDEO, BYPASS, 1.8, 2.0, 2.2, 2.4, 2.6, PACS)
5. FILTER  
Set the sharpness of image (Softest, Soft, Normal, Sharp, Sharpest)
6. OVER SCAN  
Adjust the displayed size. (0~8)
7. IMAGE SETTING  
Change the image setting. (Preset 1,2 / User 1,2,3)
8. ZOOM / PAN  
Enlarge the image, move image left and right.
9. FREEZE FRAME  
Keep the image still.

## SETUP



1. LANGUAGE  
Change the OSD language (8 language)
2. OSD COLOR  
Adjust the OSD background from white opaque to translucent.
3. OSD POSITION  
Change the OSD position. (9 Positions)
4. DURATION  
Adjust the length of time the OSD is displayed on the screen.  
(5, 10, 20, 30, 60, 90, 120, 180, 240 seconds)
5. RESET SETTING  
Change all the OSD values to factory outgoing state.
6. AUTO SOURCE SELECT  
Disable or enable auto source select.  
ON: Search through all possible input source until an active video source is found.  
OFF: Video input is manually selected.
7. USER NAME ENTRY  
Create or change the user name when powered on.


## PIP



1. LAYOUT  
Change the OSD layout. (Single, PIP, PBP1, PBP2)
2. SOURCE  
Change the secondary source.
3. SIZE  
Change the PIP size (Small, Large)
4. POSITION  
Change the secondary source.
5. SWAP  
Swap the position and size of the Primary and Secondary image.

## ▶ DVI OPTICAL / DVI DIGITAL input source

### ADJUST



#### 1. BRIGHTNESS

Increase or decrease the brightness. (Range : 0~100)


#### 2. CONTRAST

Increase or decrease the contrast. (Range : 0~100)

#### 3. BACKLIGHT

Adjust backlight dimming level. (Range : 0~100)

### COLOR TEMP



#### 1. MODE

Change the color temperature mode. C1(Reddish,6500K), C2(Bluish,9300K), USER(7200K)

#### 2. RED

Red balance (Only works with the USER mode) (Range : 0~100)


#### 3. GREEN

Green balance (Only works with the USER mode) (Range : 0~100)

#### 4. BLUE

Blue balance (Only works with the USER mode) (Range : 0~100)

### IMAGE



#### 1. IMAGE SIZE

Change the image size. (Full, Fill aspect\*\*, 1:1\*\*, Normal)

\*\* Only in DVI Optical

#### 2. GAMMA

Adjust GAMMA value (VIDEO, BYPASS, 1.8, 2.0, 2.2, 2.4, 2.6, PACS)

#### 3. FILTER

Set the sharpness of image (Softest, Soft, Normal, Sharp, Sharpest)

#### 4. OVER SCAN

Adjust the displayed size. (0~8)

#### 5. IMAGE SETTING

Change the image setting. (Preset 1,2 / User 1,2,3)


#### 6. ZOOM/PAN

Enlarge the image, move image left and right.

#### 7. FREEZE FRAME

Keep the image still.

## SETUP



### 1. LANGUAGE

Change the OSD language (8 language)

### 2. OSD COLOR

Adjust the OSD background from white opaque to translucent.

### 3. OSD POSITION

Change the OSD position. (9 Positions)

### 4. DURATION

Adjust the length of time the OSD menu is displayed on the screen.

(5, 10, 20, 30, 60, 90, 120, 180, 240 seconds)

### 5. RESET SETTING

Change all the OSD values to factory outgoing state.

### 6. AUTO SOURCE SELECT

Disable or enable auto source select.


ON: Search through all possible input source until an active video source is found.

OFF: Video input is manually selected.

### 7. USER NAME ENTRY

Create or change the user name when powered on.

## PIP



### 1. LAYOUT

Change the OSD layout. (Single, PIP, PBP1, PBP2)

### 2. SOURCE

Change the secondary source.

### 3. SIZE

Change the PIP size (Small, Large)

### 4. POSITION


Change the secondary source.

### 5. SWAP

Swap the position and size of the Primary and Secondary image.


## ▶ YPbPr input source

### ADJUST




1. BRIGHTNESS  
Increase or decrease the brightness. (Range : 0~100)
2. CONTRAST  
Increase or decrease the Contrast. (Range : 0~100)
3. SHARPNESS  
Adjust the sharpness of video image. (Range : 0~100)
4. SATURATION  
Change the tone of color. (Range : 0~100)
5. COLOR  
Change the richness of color. (Range : Greenish 0~50, Redish 0~50)
6. BACKLIGHT  
Adjust backlight dimming level. (Range : 0~100)
7. CLOCK  
Increase or decrease the sampling. (Range : 0~100)
8. PHASE  
Increase or decrease the Phase level. (Range : 0~100)

### COLOR TEMP




1. MODE  
Change the color temperature mode. C1(Reddish,6500K), C2(Bluish,9300K), USER (7200K)
2. RED  
Red balance. (Only works with USER Mode) (Range : 0~100)
3. GREEN  
Green balance. (Only works with USER Mode) (Range : 0~100)
4. BLUE  
Blue balance.(Only works with USER Mode) (Range : 0~100)

### IMAGE




1. IMAGE SIZE  
Change the image size. (Full, Fill aspect, 1:1, Normal, Anamorphic)
2. H POSITION  
Adjust the horizontal position of the displayed source image. (Range : 0~100)
3. V POSITION  
Adjust the vertical position of the displayed source image. (Range : 0~100)
4. GAMMA  
Adjust GAMMA value (VIDEO, BYPASS, 1.8, 2.0, 2.2, 2.4, 2.6, PACS)
5. FILTER  
Set the sharpness of image (Softest, Soft, Normal, Sharp, Sharpest)
6. OVER SCAN  
Adjust the displayed size. (0~8)
7. IMAGE SETTING  
Change the image setting. (Preset 1,2 / User 1,2,3)
8. ZOOM / PAN  
Enlarge the image, move image left and right.
9. FREEZE FRAME  
Keep the image still.

## SETUP



1. LANGUAGE  
Change the OSD language (8 language)
2. OSD COLOR  
Adjust the OSD background from white opaque to translucent.
3. OSD POSITION  
Change the osd position. (9 Positions)
4. DURATION  
Adjust the length of time the OSD menu is displayed on the screen.  
(5, 10, 20, 30, 60, 90, 120, 180, 240 seconds)
5. RESET SETTING  
Change all the OSD values to factory outgoing state.
6. AUTO SOURCE SELECT  
Disable or enable auto source select.  
ON: Search through all possible input source until an active video source is found.  
OFF: Video input is manually selected.
7. USER NAME ENTRY  
Create or change the user name when powered on.


## PIP



1. LAYOUT  
Change the OSD layout. (Single, PIP, PBP1, PBP2)
2. SOURCE  
Change the secondary source.
3. SIZE  
Change the PIP size (Small, Large)
4. POSITION  
Change the secondary source.
5. SWAP  
Swap the position and size of the Primary and Secondary image.

## ▶ SVIDEO / CVIDEO input source

### ADJUST



#### 1. BRIGHTNESS

Increase or decrease the brightness. (Range : 0~100)

#### 2. CONTRAST

Increase or decrease the Contrast. (Range : 0~100)

#### 3. SHARPNESS

Adjust the sharpness of video image. (Range : 0~100)

#### 4. SATURATION

Change the tone of color. (Range : 0~100)


#### 5. COLOR

Change the richness of color. (Range : Greenish 0~50, Redish 0~50)

#### 6. BACKLIGHT

Adjust backlight dimming level. (Range : 0~100)

### COLOR TEMP



#### 1. MODE

Change the color temperature mode. C1(Reddish,6500K), C2(Bluish,9300K), USER(7200K)

#### 2. RED

Red balance. (Only works with USER Mode) (Range : 0~100)


#### 3. GREEN

Green balance. (Only works with USER Mode) (Range : 0~100)

#### 4. BLUE

Blue balance. (Only works with USER Mode) (Range : 0~100)

### IMAGE



#### 1. IMAGE SIZE

Change the image size. (Full, Fill aspect, 1:1, Normal, Anamorphic)

#### 2. H POSITION

Adjust the horizontal position of the displayed source image. (Range : 0~100)

#### 3. V POSITION

Adjust the vertical position of the displayed source image. (Range : 0~100)

#### 4. GAMMA

Adjust GAMMA value (VIDEO, BYPASS, 1.8, 2.0, 2.2, 2.4, 2.6, PACS)

#### 5. FILTER

Set the sharpness of image ( Softest, Soft, Normal, Sharp, Sharpest)

#### 6. OVER SCAN

Adjust the displayed size. (0~8)

#### 7. IMAGE SETTING

Change the image setting.(Preset 1,2 / User 1,2,3)


#### 8. ZOOM / PAN

Enlarge the image, move image left and right.

#### 9. FREEZE FRAME


Keep the image still.

## SETUP



1. LANGUAGE  
Change the OSD language (8 language)
2. OSD COLOR  
Adjust the OSD background from white opaque to translucent.
3. OSD POSITION  
Change the osd position. (9 Positions)
4. DURATION  
Adjust the length of time the OSD menu is displayed on the screen.  
(5, 10, 20, 30, 60, 90, 120, 180, 240 seconds)
5. RESET SETTING  
Change all the OSD values to factory outgoing state.
6. AUTO SOURCE SELECT  
Disable or enable auto source select,  
ON: Search through all possible input source until an active video  
source is found.  
OFF: Video input is manually selected.
7. USER NAME ENTRY  
Create or change the user name when powered on.


## PIP



1. LAYOUT  
Change the OSD layout. (Single, PIP, PBP1, PBP2)
2. SOURCE  
Change the secondary source.
3. SIZE  
Change the PIP size (Small, Large)
4. POSITION  
Change the secondary source.
5. SWAP  
Swap the position and size of the Primary and Secondary image.


## SDI input source

### ADJUST




1. BRIGHTNESS  
Increase or decrease the brightness. (Range : 0~100)
2. CONTRAST  
Increase or decrease the Contrast. (Range : 0~100)
3. SHARPNESS  
Adjust the sharpness of video image. (Range : 0~100)
4. SATURATION  
Change the tone of color. (Range : 0~100)
5. COLOR  
Change the richness of color. (Range : Greenish 0~50, Redish 0~50)
6. BACKLIGHT  
Adjust backlight dimming level. (Range : 0~100)

### COLOR TEMP




1. MODE  
Change the color temperature mode. C1(Reddish,6500K), C2(Bluish,9300K), USER(7200K)
2. RED  
Red balance. (Only works with USER Mode) (Range : 0~100)
3. GREEN  
Green balance. (Only works with USER Mode) (Range : 0~100)
4. BLUE  
Blue balance. (Only works with USER Mode) (Range : 0~100)

### IMAGE




1. IMAGE SIZE  
Change the image size.(Full, Fill aspect, 1:1, Normal, Anamorphic)
2. GAMMA  
Adjust GAMMA value (VIDEO, BYPASS, 1.8, 2.0, 2.2, 2.4, 2.6, PACS)
3. FILTER  
Set the sharpness of image (Softest, Soft, Normal, Sharp, Sharpest)
4. OVER SCAN  
Adjust the displayed size. (0~8)
5. IMAGE SETTING  
Change the image setting. (Preset 1,2 / User 1,2,3)
6. ZOOM / PAN  
Enlarge the image, move image left and right.
7. FREEZE FRAME  
Keep the image still.

## SETUP



1. LANGUAGE  
Change the OSD language (8 language)
2. OSD COLOR  
Adjust the OSD background from white opaque to translucent.
3. OSD POSITION  
Change the OSD position. (9 Positions)
4. DURATION  
Adjust the length of time the OSD menu is displayed on the screen.  
(5, 10, 20, 30, 60, 90, 120, 180, 240 seconds)
5. RESET SETTING  
Change all the OSD values to factory outgoing state.
6. AUTO SOURCE SELECT  
Disable or enable auto source select.  
ON: Search through all possible input source until an active video source is found.  
OFF: Video input is manually selected.
7. USER NAME ENTRY  
Create or change the user name when powered on.













## PIP












1. LAYOUT  
Change the OSD layout. (Single, PIP, PBP1, PBP2)
2. SOURCE  
Change the secondary source.
3. SIZE  
Change the PIP size (Small, Large)
4. POSITION  
Change the secondary source.
5. SWAP  
Swap the position and size of the Primary and Secondary image.


OSD System overview

Menus	Function Descriptions
<b>BRIGHTNESS</b>	Press the BRIGHTNESS button to display the ADJUST menu or UP, DOWN hot key. Setting the brightness too high or too low will decrease the amount of visible grayscales.
<b>CONTRAST</b>	Press the CONTRAST button to display the ADJUST menu or +, - hot key. Setting the Contrast too high or too low will cause loss of some grayscales.
<b>CLOCK</b>	Do not adjust. It will adjust automatically after auto adjustment. When frequency value is wrong, the horizontal image will have a wrong size or noise.
<b>PHASE</b>	Do not adjust. It will adjust automatically after auto adjustment. When frequency value is wrong, the image will have a noise.
<b>BACKLIGHT</b>	Adjust backlight dimming level. Setting the backlight too low will cause dark image and too high will decrease the backlight life time.
<b>AUTO ADJUST</b>	Fit to the most appropriate screen on the D-SUB Analog signal
<b>SHARPNESS</b>	Adjust the sharpness of video image
<b>SATURATION</b>	Change the tone of color
<b>COLOR</b>	Change the richness of color(Range Greenish 0~50,Redish 0~50)
<b>COLOR TEMP C1</b>	Default 6500K color setting
<b>COLOR TEMP C2</b>	Default 9300K color setting
<b>COLOR TEMP USER</b>	Default 7200K color setting, but can be adjusted by the user
<b>IMAGE SIZE</b>	<p><b>DSUB/ DVI OPTICAL / DVI DIGITAL input source</b></p> <p><b>FULL</b> <b>FILL ASPECT</b></p> <p><b>1:1</b> <b>NORMAL</b></p> <p><b>CAUTION : FILL ASPECT,NORMAL</b> <b>Size depends on input size ratio</b></p> <p><b>YPbPr / RGBS / SDI / CVIDEO / SVIDEO input source</b></p> <p><b>FULL</b> <b>FILL ASPECT</b> <b>1:1</b></p> <p><b>NORMAL</b> <b>ANAMORPHIC</b></p> <p><b>CAUTION : FILL ASPECT,NORMAL</b> <b>Size depends on input size ratio</b></p>

Menus	Function Descriptions
<b>H POSITION</b>	Adjust the Horizontal position of the image. It will return to the default state when executing AUTO ADJUST OR RESET SETTINGS. 
<b>V POSITION</b>	Adjust the Vertical position of the image. It will return to the default state when executing AUTO ADJUST OR RESET SETTINGS. 
<b>GAMMA</b>	Adjust the gamma curve of video image.  (NOTICE) BYPASS is depends on panel gamma value, please refer to the panel specification.
<b>FILTER</b>	Adjust the gamma curve of video image.  
<b>OVER SCAN</b>	Enable 10% over scan of original input image. 0:  1:  2:  3:  4:  5:  6:  7:  8: 

Menus	Function Descriptions
<b>IMAGE SETTING</b>	Save 5 user setting values for BRIGHTNESS,CONTRAST,COLOR,TEMP,FILTER separately.(PRESET1,2,3 / USER 1,2)
<b>ZOOM / PAN</b>	<p>Controls the zoom in/out of the image.</p> <p>0 Step :  4 Step :  9 Step : </p> <p>13 Step :  18 Step :  22 Step : </p> <p>26 Step :  30 Step : </p> <p>Controls the PAN in/out of the image.</p> 
<b>FREEZE FRAME</b>	Freeze the main image. Does not support secondary image at PIP mode.
<b>LANGUAGE</b>	Change the OSD language (8 languages) ENGLISH / GERMAN / FRENCH / SPANISH / ITALIAN / JAPANESE / CHINESE / KOREAN
<b>OSD COLOR</b>	Adjust the color of the OSD.
	
<b>OSD POSITION</b>	Adjust the position of the OSD.
	

Menus	Function Descriptions																																																																																									
<b>OSD DURATION</b>	Adjust the length of time the OSD menu is displayed on the screen.																																																																																									
<b>RESET SETTING</b>	Reset the unit to factory outgoing status.																																																																																									
<b>AUTO SOURCE SELECT</b>	When AUTO SOURCE SELECT is on, the monitor will automatically search input source except PIP sub Source.																																																																																									
<b>USER NAME ENTRY</b>	Enter user name or other information.																																																																																									
<b>PIP LAY OUT</b>	Change sub window layout (SINGLE, PIP, PBP1, PBP2) PIP : Main/Sub window doesn't change aspect ratio PBP1 : Main/Sub window doesn't change aspect ratio PBP2 : Main/Sub window H/V display 1:1 size																																																																																									
<b>PIP SOURCE</b>	Select PIP source input. You can change other sub windows through the OSD PIP menu. Below chart is PIP																																																																																									
	<table border="1"> <thead> <tr> <th rowspan="2">Input Source</th> <th colspan="8">Sub window</th> </tr> <tr> <th>DVI OPTICAL</th> <th>DVI DIGITAL</th> <th>DSUB ANALOG</th> <th>SDI</th> <th>YPbPr</th> <th>RGBS</th> <th>CVIDEO</th> <th>SVIDEO</th> </tr> </thead> <tbody> <tr> <td>DVI OPTICAL</td> <td>X</td> <td>O</td> <td>O</td> <td>O</td> <td>O</td> <td>O</td> <td>O</td> <td>O</td> </tr> <tr> <td>DVI DIGITAL</td> <td>O</td> <td>X</td> <td>X</td> <td>O</td> <td>X</td> <td>X</td> <td>O</td> <td>O</td> </tr> <tr> <td>DSUB ANALOG</td> <td>O</td> <td>X</td> <td>X</td> <td>O<sup>1</sup></td> <td>X</td> <td>X</td> <td>O</td> <td>O</td> </tr> <tr> <td>SDI</td> <td>O</td> <td>O</td> <td>O<sup>1</sup></td> <td>X</td> <td>O<sup>1</sup></td> <td>O<sup>1</sup></td> <td>X</td> <td>X</td> </tr> <tr> <td>YPbPr</td> <td>O</td> <td>X</td> <td>X</td> <td>O<sup>1</sup></td> <td>X</td> <td>X</td> <td>O</td> <td>O</td> </tr> <tr> <td>RGBS</td> <td>O</td> <td>X</td> <td>X</td> <td>O<sup>1</sup></td> <td>X</td> <td>X</td> <td>O</td> <td>O</td> </tr> <tr> <td>CVIDEO</td> <td>O</td> <td>O</td> <td>O</td> <td>X</td> <td>O</td> <td>O</td> <td>X</td> <td>X</td> </tr> <tr> <td>SVIDEO</td> <td>O</td> <td>O</td> <td>O</td> <td>X</td> <td>O</td> <td>O</td> <td>X</td> <td>X</td> </tr> </tbody> </table> <p>O<sup>1</sup>: Support up to UXGA,60Hz (162MHz)</p>	Input Source	Sub window								DVI OPTICAL	DVI DIGITAL	DSUB ANALOG	SDI	YPbPr	RGBS	CVIDEO	SVIDEO	DVI OPTICAL	X	O	O	O	O	O	O	O	DVI DIGITAL	O	X	X	O	X	X	O	O	DSUB ANALOG	O	X	X	O <sup>1</sup>	X	X	O	O	SDI	O	O	O <sup>1</sup>	X	O <sup>1</sup>	O <sup>1</sup>	X	X	YPbPr	O	X	X	O <sup>1</sup>	X	X	O	O	RGBS	O	X	X	O <sup>1</sup>	X	X	O	O	CVIDEO	O	O	O	X	O	O	X	X	SVIDEO	O	O	O	X	O	O	X	X
Input Source	Sub window																																																																																									
	DVI OPTICAL	DVI DIGITAL	DSUB ANALOG	SDI	YPbPr	RGBS	CVIDEO	SVIDEO																																																																																		
DVI OPTICAL	X	O	O	O	O	O	O	O																																																																																		
DVI DIGITAL	O	X	X	O	X	X	O	O																																																																																		
DSUB ANALOG	O	X	X	O <sup>1</sup>	X	X	O	O																																																																																		
SDI	O	O	O <sup>1</sup>	X	O <sup>1</sup>	O <sup>1</sup>	X	X																																																																																		
YPbPr	O	X	X	O <sup>1</sup>	X	X	O	O																																																																																		
RGBS	O	X	X	O <sup>1</sup>	X	X	O	O																																																																																		
CVIDEO	O	O	O	X	O	O	X	X																																																																																		
SVIDEO	O	O	O	X	O	O	X	X																																																																																		
<b>PIP SIZE</b>	Change PIP mode sub window size  <div style="display: flex; justify-content: space-around;"> <span>SMALL</span> <span>LARGE</span> </div> 																																																																																									
<b>PIP POSITION</b>	Change PIP mode sub window position  																																																																																									

Menus	Function Descriptions
<b>SWAP</b>	<p>Change main window and sub window position on PIP,PBP1,PBP2.</p>  <p>The diagram illustrates three window configurations for the SWAP function:</p> <ul style="list-style-type: none"> <li><b>PIP:</b> Shows a large "Main" window at the bottom and a smaller "Sub" window at the top right. A double-headed arrow between these two states is labeled "PIP".</li> <li><b>PBP1:</b> Shows a horizontal bar divided into two equal sections, with "Main" on the left and "Sub" on the right. A double-headed arrow between these two states is labeled "PBP1".</li> <li><b>PBP2:</b> Shows a horizontal bar divided into two equal sections, with "Sub" on the left and "Main" on the right. A double-headed arrow between these two states is labeled "PBP2".</li> </ul>

## Standard Signal table

### ▶ PC Supported Mode

Resolution	Horizontal Frequency (KHz)	Vertical Frequency (Hz)	Clock Frequency (MHz)			
640 X 350 @70Hz	31.50	70.00	25.17			
720 X 400 @70Hz	31.50	70.00	28.32			
640 X 480 @60Hz	31.50	60.00	25.18			
640 X 480 @75Hz	37.50	75.00	31.50			
800 X 600 @60Hz	37.90	60.00	40.00			
800 X 600 @75Hz	46.90	75.00	49.50			
1024 X 768 @60Hz	48.40	60.00	65.00			
1024 X 768 @75Hz	60.00	75.00	78.75			
1152 X 864 @60Hz	54.35	60.00	80.000			
1152 X 864 @75Hz	67.50	75.00	108.00			
1280 X 720 @60Hz	44.70	60.00	74.40			
1280 X 1024@60Hz	64.00	60.00	108.00			
1280 X 1024@75Hz	80.00	75.00	135.00			
1360 X 768@75Hz	47.70	60.00	108.75			
1600 X 1200@60Hz	75.00	60.00	162.00			
1920 X 1080@60Hz	67.50	60.00	148.50			
1920 X 1200@60Hz	74.00	60.00	154.12			


### ▶ SDI Video format

Output Signal	Description
SMPTE-274M	1080i (60 / 59.94 / 50) 1080p (30 / 29.97 / 25 / 24 / 24sF / 23.98 / 23.98sF)
SMPTE-296M	720p (60 / 59.94 / 50)
SMPTE-260M	1035i (60 / 59.94)
SMPTE-125M	480i (59.94)
ITU-R BT.656	576i (50)

## Signal connector Pin Assignments


### VGA (15Pin D-SUB)

Pin No.	Assignment	Pin No.	Assignment
1	Red	9	No Connection
2	Green	10	Ground-Sync
3	Blue	11	Ground
4	Ground	12	DDC Data
5	DDC 5V Standby Cable Connection check	13	H.Sync
		14	V.Sync
6	Ground-Red	15	DDC Clock
7	Ground-Green		
8	Ground-Blue		




### DVI In,Out (24Pin DVI-D)

Pin No.	Assignment	Pin No.	Assignment
1	T.M.D.S. Data2-	13	No Connection
2	T.M.D.S. Data2+	14	+5V Power
3	T.M.D.S. Data2 Shield	15	Ground
4	No Connection	16	Hot Plug Detect
5	No Connection	17	T.M.D.S. Data0-
6	DDC Clock	18	T.M.D.S. Data0+
7	DDC Data	19	T.M.D.S. Data0 Shield
8	No Connection	20	No Connection
9	T.M.D.S. Data1-	21	No Connection
10	T.M.D.S. Data1+	22	T.M.D.S. Clock Shield
11	T.M.D.S. Data1 Shield	23	T.M.D.S. Clock+
12	No Connection	24	T.M.D.S. Clock-




### C-Video (BNC)

Pin No.	Assignment
1	Composite
2	Ground




▶ S-Video (BNC)

Pin No.	Assignment
1	S-VIDEO/Y (Luma)
2	S-VIDEO/C (Chroma)
3	Ground




▶ RS232C (D-SUB 9Pin)

Pin No.	Assignment
1	No Connection
2	TXD
3	RXD
4	No Connection
5	Ground
6	No Connection
7	No Connection
8	No Connection
9	No Connection




▶ SDI (BNC)

Pin No.	Assignment
1	SDI IN
2	SDI OUT
3	Ground




### ▶ RGBHV/RGBS/YPbPr (BNC)

Pin No.	Assignment	
	RGBS	Y Pb Pr
1	Red	Pr
2	Green	Y
3	Blue	Pb
4	H-Sync / C-Sync	No Connection
5	V-Sync	No Connection
6	Ground	




### ▶ OPTICAL

Pin No.	Assignment
1	OPTICAL Clock
2	OPTICAL Blue
3	OPTICAL Green
4	OPTICAL Red



### ▶ GPIO

Pin No.	Assignment
1	P,S Swap
2	PIP,PBP1,PBP2 Select
3	Record Indicator
4	Ground



## Specification

▶ FS-L190\*D / FS-L190\*DT

Model		<b>FS-L190*D / FS-L190*DT</b>
Optical Characteristics	Type	19" TFT-LCD
	Screen Size	19 inch
	Maximum Resolution	1280 X 1024 @ 60Hz
	Pixel Pitch	0.294(H) mm X 0.294(V) mm
	Display Colors	16.7M
	Contrast Ratio(Typical)	500:1
	Viewing Angle	85° / 85° / 85° / 85°
	Response Time	14 msec(Gray to Gray)
	Luminance(Typical)	210cd/m <sup>2</sup>
Touch Specification (F-L190*DT only)	Touch panel	ELO 5wire resistive touch screen
	Interface	USB & SERIAL
Front Filter (FS-L190*D only)	Acrylic	Double side Anti-Reflection coating
Synchronization	Horizontal Frequency	30KHz~75KHz
	Vertical Frequency	50Hz~75Hz
Power Consumption	Maximum	Max 60W
	Standby mode	Max 20W
Control Key	Front side	INPUT, -, +, ▼, ▲, PIP, MENU, POWER
Input Signal	Video	1개 DVI, 1개 Optical DVI 1(DVI 2 option), 1개 D-SUB, 1개 BNC (CVBS)Video, 2 개 BNC (SVHS Y/C), 1 개 BNC(SDI), 5 개 BNC (Component Y/G, Pb/B, Pr/R, H/CS, VS)
Output Signal	Video	1개 DVI, 1개 BNC (SDI)
Input Power	DC 12V, 7A Max	
Dimension	Size and Weight	423(W)X351.5(H)X76.5(D) (mm) 7.3Kg -Without stand 16.653(W)X13.838(H)X3.011(D) (inch) 16.09 lbs-Without stand

▶ FS-L240\*D / FS-L240\*DT

Model	<b>FS-L240*D / FS-L240*DT</b>	
Optical Characteristics	Type	24" TFT-LCD
	Screen Size	24 inch
	Maximum Resolution	1920 X 1200 @ 60Hz
	Pixel Pitch	0.270(H) mm X 0.270(V) mm
	Display Colors	16.7M
	Contrast Ratio(Typical)	700:1
	Viewing Angle	89° / 89° / 89° / 89°
	Response Time	5 msec(Gray to Gray)
	Luminance(Typical)	320cd/m <sup>2</sup>
Touch Specification (F-L240*DT only)	Touch Panel	ELO 5wire resistive touch screen
	Interface	USB & SERIAL
Front Filter (FS-L240*D only)	Acrylic	Double side Anti-Reflection coating
Synchronization	Horizontal Frequency	30KHz~75KHz
	Vertical Frequency	50Hz~75Hz
Power Consumption	Maximum	Max 100W
	Standby Mode	Max 20W
Control Key	Front side	INPUT, -, +, ▼, ▲, PIP, MENU, POWER
Input Signal	Video	1개 DVI, 1개 Optical DVI 1(DVI 2 option), 1개 D-SUB, 1개 BNC (CVBS)Video, 2 개 BNC (SVHS Y/C), 1 개 BNC(SDI), 5 개 BNC (Component Y/G, Pb/B, Pr/R, H/CS, VS)
Output Signal	Video	1개 DVI, 1개 BNC (SDI)
Input Power	DC 24V, 6.25A Max	
Dimension	Size and Weight	580(W)X386(H)X95(D) (mm) 7.5Kg -Without stand 22.834(W)X15.196(H)X3.740(D) (inch) 16.53 lbs-Without stand

▶ FS-L260\*D

Model		<b>FS-L260*D</b>
Optical Characteristics	Type	26" TFT-LCD
	Screen Size	26 inch
	Maximum Resolution	1920X 1200 @ 60Hz
	Pixel Pitch	0.2685(H) mm X 0.2685(V) mm
	Display Colors	16.7M
	Contrast Ratio(Typical)	700:1
	Viewing Angle	89° / 89° / 89° / 89°
	Response Time	5 msec(Gray to Gray)
	Luminance(Typical)	320cd/ m <sup>2</sup>
Synchronization	Horizontal Frequency	30KHz~75KHz
	Vertical Frequency	50Hz~75Hz
Front Filter	Acrylic	Double side Anti-Reflection coating
Power Consumption	Maximum	Max 130W
	Standby Mode	Max 20W
Control Key	Front side	INPUT, -, +, ▼, ▲, PIP, MENU, POWER
Input Signal	Video	1개 DVI, 1개 Optical DVI 1(DVI 2 option), 1개 D-SUB, 1개 BNC (CVBS)Video, 2 개 BNC (SVHS Y/C), 1 개 BNC(SDI), 5 개 BNC (Component Y/G, Pb/B, Pr/R, H/CS, VS)
Output Signal	Video	1개 DVI, 1개 BNC (SDI)
Input Power	DC 24V, 6.25A Max	
Dimension	Size and Weight	618(W)X412(H)X99.5(D) (mm) 8.8Kg -Without stand 24.330(W)X16.220(H)X3.917(D) (inch) 19.4bs-Without stand

▶ FS-L320\*D

Model	FS-L320*D	
Optical Characteristics	Type	32" TFT-LCD
	Screen Size	32 inch
	Maximum Resolution	1920X 1080 @ 60Hz
	Pixel Pitch	0.363(H) mm X 0.363(V) mm
	Display Colors	1.06 billion
	Contrast Ratio(Typical)	850:1
	Viewing Angle	89° / 89° / 89° / 89°
	Response Time	6msec(Gray to Gray)
	Luminance(Typical)	380cd/m²
Synchronization	Horizontal Frequency	30KHz~75KHz
	Vertical Frequency	50Hz~75Hz
Front Filter	Glass	Double side Anti-Reflection coating
Power Consumption	Maximum	Max 150W
	Standby Mode	Max 20W
Control Key	Front Side	INPUT, -, +, ▼, ▲, PIP, MENU, POWER
Input Signal	Video	1개 DVI, 1개 Optical DVI 1(DVI 2 option), 1개 D-SUB, 1개 BNC (CVBS)Video, 2 개 BNC (SVHS Y/C), 1 개 BNC(SDI), 5 개 BNC (Component Y/G, Pb/B, Pr/R, H/CS, VS)
Output Signal	Video	1개 DVI, 1개 BNC (SDI)
Input Power	DC 24V, 7.5A Max	
Dimension	Size and Weight	795(W)X485(H)X100(D) (mm) 15.8Kg -Without stand 31.299(W)X19.094(H)X3.937(D) (inch) 34.8lbs-Without stand

▶ FS-L420\*D

Model	<b>FS-L420*D</b>	
Optical Characteristics	Type	42" TFT-LCD
	Screen Size	42 inch
	Maximum Resolution	1920X 1080 @ 60Hz
	Pixel Pitch	0.4845(H) mm X 0.4845(V) mm
	Display Colors	1.06 billion
	Contrast Ratio(Typical)	650:1
	Viewing Angle	89° / 89° / 89° / 89°
	Response Time	5 msec(Gray to Gray)
	Luminance(Typical)	400cd/ m <sup>2</sup>
Synchronization	Horizontal Frequency	30KHz~75KHz
	Vertical Frequency	50Hz~75Hz
Front Filter	Glass	Double side Anti-Reflection coating
Power Consumption	Maximum	Max 260W
	Standby Mode	Max 20W
Control Key	Front Side	INPUT, -, +, ▼, ▲, PIP, MENU, POWER
Input Signal	Video	1 個 DVI, 1 個 Optical DVI 1(DVI 2 option), 1 個 D-SUB, 1 個 BNC (CVBS)Video, 2 個 BNC (SVHS Y/C), 1 個 BNC(SDI), 5 個 BNC (Component Y/G, Pb/B, Pr/R, H/CS, VS)
Output Signal	Video	1 個 DVI, 1 個 BNC (SDI)
Input Power	AC 100~230V, 50~60Hz, 3A Max	
Dimension	Size and Weight	1024.6(W)X617.4(H)X111.1(D) (mm) 28Kg -Without stand 40.338(W)X24.307(H)X4.374(D) (inch) 61.7 lbs-Without stand